

Soldier 2.0




High Flow Microcatheter with Ultra-Thin Wall Technology




Largest Inner
Diameter
0.022" ID
2.0 Fr OD

Fastest Flow
Rate at 1000 PSI




Best Embolic
Compatibility




Exceptional
Tracking



High
Torquability



Unparalleled
Performance




Best in Class Microcatheter

Mastering
Performance
and Delivery

Soldier microcatheters are designed with advanced technology, allowing them to have the highest level of performance of any microcatheter available. Tracking, torquing, curved tips, largest lumen, highest flow rate, high-pressure tolerance...



...get where you want **FAST**, deliver what you want **ALWAYS**.



Unique 2.0 Fr profile with 0.022" inner diameter
Largest ID, Highest Flow Rate, Best Embolic Compatibility

Product Code (part number)	Outer Diameter (Distal)	Outer Diameter (Proximal)	Working Length	Inner Diameter	Tip Shape
S22-SP-130	2.0 Fr	2.0 Fr	130 cm	0.022"	Straight
S22-SP-155	2.0 Fr	2.0 Fr	155 cm	0.022"	Straight
S22-KP-130	2.0 Fr	2.0 Fr	130 cm	0.022"	K-Tip
S22-KP-155	2.0 Fr	2.0 Fr	155 cm	0.022"	K-Tip



INDICATIONS FOR USE

The Soldier is indicated for use in the blood vessels of the peripheral vasculature. It is intended to assist in the delivery of diagnostic agents and therapeutic agents into the target treatment area.

CONTRAINDICATIONS

Not intended for use in embolectomy or angioplasty procedures. Not intended for use in the coronary vessels. Not intended for use in the neurovasculature.

Embolx, Inc.
530 Lakeside Dr. #150
Sunnyvale, CA 94085
www.embolx.com

To order product or for more information
contact us at customer care@embolx.com
or +1-408-990-2949.

©2025 Embolx, Inc. All rights reserved.
Soldier is a registered trademark of Embolx, Inc.